




MC2 Viewfinder Trainer

Practice aiming the MC2 emitter with your phone.



How to use

1. Lay this page flat on a table.
2. Open the trainer (scan the QR on the right).
3. Allow camera access when prompted.
4. Point your phone at the marker like the MC2 emitter.
5. Tap Tutorial to learn, or Play to free-practice.



Open the trainer
(scan with your phone)

ArUco DICT_4X4_1000 ID 0 - marker 6 cm x 6 cm

Built by the Advanced Concepts team - OXOS Medical